

The **Real World MDT** multiple drug test kit is a rapid, one-step lateral flow immunoassay in cassette format for the qualitative detection of AMP, COC, THC metabolites in urine, with cut-off levels for the most specific and sensitive detection of drug usage.

### Features

One Step  
 Highly Accurate  
 Clear Results  
 Results after 5 Minutes  
 Test anytime of the day



### Test Performance

Quality Specifications:

- ✓ Accuracy- 98 – 99%
- ✓ Sensitivity- 95 – 100%
- ✓ Specificity- 99 – 100%
- ✓ Reproducibility- 99%+

The Real World MDT is manufactured under ISO 13485:2003 quality assurance conditions, and is a single-step immunoassay tests for the qualitative detection of drugs in urine.

### Summary and Test Cut-Off levels

Can detect drug metabolites in urine up to 72 hours after drug usage, at the following cut-off levels:

- AMP (d-Amphetamine) 1000 ng/ml → *Tik (crystal meth, glass, speed, ice)*
- COC (Benzoylecgonine) 300 ng/ml → *Cocaine (crack, coke, rocks, snow)*
- THC (Marijuana) 50 ng/ml → *Marijuana (dagga, weed, pot, dope, hashish)*

The device utilizes colloidal gold as the detection system in a competitive binding process, and is simple to run, easy to interpret and deliver quick, efficiently and effectively.

### Instructions for Use and Interpretation of Results

Detailed instructions and interpretation of results are provided in the product insert, and as part of the product packaging

### Contact Us:

Real World Diagnostics PO BOX 38140 Howard Place, Pinelands, 7430  
 ☎ (021) 442 5385 📠 (021) 447 6219 ✉ info@real-world.co.za